

SURGICAL STERILE GOWN XREINFORCED

MZH-SC054

- These surgical gowns have the ability to resist high fluid exposure and blood-borne pathogens.
- 54 GSM
- Sterile
- SMS/PE
- Back Ties
- Neck Ties



| | |
|---|--|
| ANSI/AAMI PB70 BARRIER PERFORMANCE | AAMI Level 4 |
| TEST METHOD | ASTM F1671, Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System |
| TEST DEFINITION | ASTM F1671 Measures the resistance of materials used in protective clothing to penetration by blood-borne pathogens using a surrogate microbe under conditions of continuous liquid contact. |
| ANTICIPATED FLUID EXPOSURE | High fluid levels |

| LEVEL GUIDE | | | |
|--------------------|----------------|---------------------|-----------------|
| 1 | 2 | 3 | 4 |
| Minimal Protection | Low Protection | Moderate Protection | High Protection |

Contact your sales rep for more information or to place an order.